




EASY, FAST, RELIABLE

ELDONCARD® 2521 - for ABO & RhD forward blood typing



Recommended for:

- Medical Market – Blood Banks, Hospitals, Clinics etc.
- Emergency & Rescue – NGO's
- Military, Police and Civil Service
- Donor Recruitment
- Rhesus tests for Pregnant Women
- Educational Market
- Health / Home Market
- Bed side test

ELDONCARD® 2521 carries 2 x 4 circles with a diameter of 20 mm.

The card is used for the detection of the blood type of one/two person/s
Control field to detect false possible false positives.

The procedure is described in detail in the "Instruction for Use".

Instructional video: www.eldoncard.com (English/Espanol/French + prof. use)

Available as:



- **ELDON BAG 2521**

Contains:

- 25 ELDONCARD 2521
- 25 Blood Group Determinations

- **ELDON BOX 2521**

Contains:*

- 100 ELDONCARDS
- 2521
- 100 Blood Group Determinations



- **Costum made Kits for special use NGO/Doctors/ Military**

**Each EldonCard is in an aluminum envelope.*

Labelling: CE marked & FDA cleared and in other countries

Storage: 5-37 Celsius / 41-99 Fahrenheit – **NO COOL**

CHAIN Expiration: 24 month from first day of production

Accessories:

In general, accessories are available in numbers matching the number of Blood Group Determinations.



- Safety Lancet



- Steel Lancet



- EldonStick



- Alcohol Swabs



- EldonFoil



- Pippette

General: To perform one Blood Group Determination 2521 you need 1+1 lancet, 4x2 EldonSticks, 1 EldonFoil, 1 pipette and 1+1 alcohol swab.



ELDON
BIOLOGICALS A/S